

GPDT-GENERAL PURPOSE/Distal Traction

- The liner is generally known as the GPDT liner (General Purpose Distal Traction)
- It is manufactured in the original soft and tacky gel
- Indicated for the elderly or low activity patient seeking a soft cushioning gel with good skin adherence
- Tackiness remains throughout the lifetime of the liner
- Note: It may be used with a high activity patient BUT then two GPDT liners should be used per 18 month period. Each liner should be rested for a day after a day of use, alternate back to front as well.
- Patient instruction: when white spot develop seek O&P attention to relieve the pressure area.

