

ENCP16 / SKIN RELIEVER-GEL SHEATH

- The ALPS Skin Reliever eliminates shear forces, friction, and prevents abrasions on the skin of the stump. It offers both comfort and durability.
- It is used to accommodate residuum volume fluctuations in both elevated vacuum systems or locking liners. It is applied directly onto the skin before the application of the locking liner.
- Superior elongation allows for very little shear forces or compression against the skin.
- ALPS Skin Reliever can also be worn under silicone liners and may be used with standard P-lite liners.

